


Anterior


Posterior

ANTERIOR CROWN PREPARATION




UNIFORM FACIAL PREPARATION

A medium grit, round-ended, diamond bur is used to remove a uniform thickness of facial enamel by joining the depth-cut grooves.



INCISAL PREPARATION

The diamond bur is angled to bevel back the incisal edge.




CHAMFER MARGINS

Correct preparation of the chamfer margins interproximally allows the appropriate bulk of porcelain.




3-UNIT BRIDGE PREPARATION




100% CUSTOMER SATISFACTION
GUARANTEED!


POSTERIOR CROWN PREPARATION



CONVENTIONAL CEMENTATION PREPARATION



3-UNIT BRIDGE PREPARATION



When layered or pressed ceramic margins are preferred in conjunction with a zirconia framework, enhanced gingival esthetics can be achieved with a definitive 90 degree shoulder preparation.